




## A-010

### Bulk Enzyme Production

Arylacylamidase catalyzes the conversion of p-nitroacetanilide to p-nitroaniline, which absorbs strongly at 405 nm:



Arylacylamidase also catalyzes the conversion of acetaminophen to p-aminophenol:



In one method, p-aminophenol reacts with o-cresol in the presence of periodate to form an indophenol, which absorbs strongly at 615 nm.

Used for the enzymatic determination of acetaminophen in serum

Measures the increase of p-aminophenol via secondary colorimetric chemistries

# A-010 Bulk Enzyme Production

## Specifications

### Form

White, off-white, or yellow lyophilized powder.

### Activity

≥15 U/mg powder.

### Unit

One unit is defined as the amount of enzyme which catalyzes the conversion of 1 μmole of p-nitroacetanilide to p-nitroaniline at 30°C, pH 8.5 under the conditions given in the assay procedure.

## Assay Method

### Reagents

- 1 Tris-HCl buffer: 50 mM, pH 7.0.
- 2 Stock Substrate Solution: Dissolve 3.03 g of Tris Base in 450 mL of DI water. Heat to 65 ± 5°C. With stirring, add 90.1 mg of p-nitroacetanilide and stir vigorously until it is dissolved. Do not over-heat. Cool the solution to 24 ± 2°C but not less than 20°C. Adjust solution to pH 8.50 ± 0.02 with 5M HCl. Absorbance at 405 nm must be ≤0.2 vs. DI water. Store in amber bottle at room temperature.
- 3 Diluted β-Mercaptoethanol: Add 50 μL of β-mercaptoethanol to 450 μL of DI water. Prepare fresh daily and store capped at room temperature.
- 4 Enzyme Diluent: Add 140 μL of Diluted β-mercaptoethanol to 100 mL of 50 mM Tris-HCl buffer, pH 7.0. Make fresh daily and store cold in a sealed container.
- 5 Working Reagent: Add 140 μL of Diluted β-mercaptoethanol to 100 mL of Stock Substrate Solution. Prepare fresh daily and store capped at room temperature in an amber bottle.
- 6 Enzyme Solution: Prepare a 10 mg/mL enzyme solution in enzyme diluent. Dilute the enzyme in same to yield an activity of approximately 0.10 to 0.15 U/mL.

### Procedure

Combine 0.48 mL of Working Reagent equilibrated to 30°C with 20 μL of diluted enzyme solution in a cuvette.

Mix and measure the rate of increase in absorbance between one and three minutes at 405 nm in a spectrophotometer controlled at 30°C.

The change in absorbance should be between 0.03 and 0.07 per minute.

## Properties

### Solubility

Arylacylamidase is soluble in water and buffers.

### Effect of Buffers


Arylacylamidase is stable in pH 8.5 Tris and Phosphate buffers with a molarity of 20 mM to 100mM. The following chart shows the percent activity obtained when assayed in various Tris and Phosphate buffers.

Concentration	Tris	Phosphate
20 mM, pH 8.5	98%	99%
50 mM, pH 8.5	100%	97%
100 mM, pH 8.5	94%	94%


### Optimum pH and Temperature

The graphs below show the relative activity of arylacylamidase at various temperatures and pH under the assay conditions with the p-nitroacetanilide substrate:

A-010 Temperature Effect



A-010 pH Effect



### Michaelis-Menten Constant

Arylacylamidase has an apparent  $K_M$  of ≥200 mM for acetaminophen and 38 μM for p-nitroacetanilide.

### pI


Arylacylamidase has an apparent pI of 4.9-5.2.

### Molecular Weight

The molecular weight of arylacylamidase was determined to be ≈ 52 kDa via size exclusion chromatography and subsequent enzyme analysis.

The image below demonstrates the electrophoretic separation of a sample from a lot of arylacylamidase. Protein standard markers are shown on the left.

The major protein migrates as a single polypeptide chain of 54 kDa. Mass spectrophotometric analysis confirms the presence of a 51 kDa protein.



### Calculation

Calculate arylacylamidase activity as follows:

$$\frac{\text{U/mg} = \Delta A_{405} \times \text{cv} \times \text{dilution}}{10 \times \text{sv}}$$

where,  
cv = reaction volume in mL  
sv = enzyme sample volume in mL

## EKF Life Sciences

1814 Leer Drive | 6879 Enterprise Drive, Suite 500  
Elkhart, IN 46514 | South Bend, IN 46628

Phone: (574) 264-7384

Toll Free: (800) 545-4437

Fax: (574) 266-0062

Email: lifesciences@ekfdiagnostics.com

