



## STAT-Site WB

Reliable and efficient  
 $\beta$ -ketone and glucose  
measurement from  
whole blood

- ◆ CLIA waived
- ◆ Results in just 5-10 seconds
- ◆ Handheld
- ◆ Easy to use

# Simple, accurate and efficient $\beta$ -ketone and glucose measurement

Simple and intuitive	<ul style="list-style-type: none"> <li>• Auto start when the strip is inserted</li> <li>• Auto reaction time countdown</li> <li>• Auto switch-off after 3 minutes without action</li> <li>• Touch button to eject test strip</li> <li>• Two 1.5V AAA alkaline batteries</li> <li>• 400 results stored on the system</li> </ul>
Accurate with excellent linear range	<ul style="list-style-type: none"> <li>• <b>Measurement range:</b>  <math>\beta</math>-ketone: 0.1 to 8.0 mmol/L                      Blood glucose: 10 to 600 mg/dL (0.5 to 33.3 mmol/L)</li> <li>• <b>Hematocrit range:</b>  <math>\beta</math>-ketone: 10 to 70%                      Blood glucose: 20 to 70%</li> <li>• Calibrated using the gold standard Liquicolor beta-Hydroxybutyrate method</li> </ul>
Fast and efficient sample handling	<ul style="list-style-type: none"> <li>• 1 <math>\mu</math>L sample size capillary or venous whole blood</li> <li>• <math>\beta</math>-ketone result in 10 seconds</li> <li>• Glucose result in 5 seconds</li> </ul>
Ergonomic design	<ul style="list-style-type: none"> <li>• Compact: 95.5 (l) x 57 (w) x 15 (h) mm, 58 g</li> <li>• Simple, large LCD screen interface</li> </ul>

Catalog No.	Description
7000-001	STAT-Site WB Dual Analyte Measurement System
7010-025	Blood Glucose Strip box 25pcs
7010-050	Blood Glucose Strip box 50pcs
7010-100	Blood Glucose Strip box foil pack 100pcs
7020-025	$\beta$ -Ketone Strip box 25pcs
7020-050	$\beta$ -Ketone box 50pcs
7020-100	$\beta$ -Ketone box foil pack 100pcs
7014-204	Blood Glucose Control Solution 1 x W2 (norm), 1 x B3 (high)
7024-204	$\beta$ -Ketone Control solution 1 x L1, 1 x L2



1 Insert test strip into the meter



2 Add blood sample



3 Result within 10 seconds for  $\beta$ -ketone and 5 seconds for glucose

Distributed by

EKF USA  
 1261 North Main Street  
 Boerne, Texas 78006  
 ☎ 830.249.0772  
 Tel (USA Toll Free): 1.800.531.5535

[www.ekfusa.com](http://www.ekfusa.com)  
[sales2@ekfdiagnostics.com](mailto:sales2@ekfdiagnostics.com)

**EKF** | Diagnostics  
 for life